

产品名称: PEX11B Rabbit Monoclonal Antibody
产品货号: AMRe84820



产品概述 (Summary)

产品名称 (Production Name)	PEX11B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PEX11B
别名 (Alternative Names)	Peroxin 11B; PEX11 beta; PEX11B; Pex11p beta; PEX14B;; Peroxin 11B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O96011.A synthesized peptide derived from human Peroxin 11B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa ; Observed MW: 25 kDa

研究背景 (Background)

Involved in peroxisomal proliferation. May regulate peroxisomes division by recruiting the dynamin-related GTPase DNM1L


产品名称: PEX11B Rabbit Monoclonal Antibody
产品货号: AMRe84820



to the peroxisomal membrane.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PEX11B expression in HT-29 cell lysate.

注意事项 (Note)

For research use only .