

产品名称: Nup107 Rabbit Monoclonal Antibody
产品货号: AMRe84821



产品概述 (Summary)

产品名称 (Production Name)	Nup107 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nup107
别名 (Alternative Names)	NUP107; NUP84;;NUP107
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P57740.A synthesized peptide derived from human NUP107

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	106 kDa

研究背景 (Background)

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Required for the assembly of peripheral


产品名称: Nup107 Rabbit Monoclonal Antibody
产品货号: AMRe84821



proteins into the NPC. May anchor NUP62 to the NPC.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Nup107 expression in 293 cell lysate.

注意事项 (Note)

For research use only .