

产品名称: PREB Rabbit Monoclonal Antibody
产品货号: AMRe84822



产品概述 (Summary)

产品名称 (Production Name)	PREB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PREB
别名 (Alternative Names)	Preb; SEC12;;PREB
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9HCU5.A synthesized peptide derived from human PREB

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	45 kDa

研究背景 (Background)

Guanine nucleotide exchange factor (GEF) that regulates the assembly of the coat protein complex II/COPII in endoplasmic


产品名称: PREB Rabbit Monoclonal Antibody
产品货号: AMRe84822



reticulum (ER) to Golgi vesicle-mediated transport. Selectively activates SAR1A and SAR1B by promoting the exchange of guanosine diphosphate (GDP) for guanosine triphosphate (GTP) in these small GTPases.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PREB expression in Raji cell lysate.

注意事项 (Note)

For research use only .