

产品名称: Cleaved-Caspase 3 p12 Rabbit Monoclonal Antibody
产品货号: AMRe84828



产品概述 (Summary)

产品名称 (Production Name)	Cleaved-Caspase 3 p12 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cleaved-Caspase 3 p12
别名 (Alternative Names)	CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1
基因 ID (Gene ID)	836.0
蛋白 ID (SwissProt ID)	P42574.Recombinant protein of human Caspase-3 p12

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 32 kDa; Observed MW: 32, 12 kDa

产品名称: Cleaved-Caspase 3 p12 Rabbit Monoclonal Antibody
产品货号: AMRe84828




研究背景 (Background)

Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce 2 subunits, large and small, that dimerize to form the active enzyme.

研究领域 (Research Area)

Apoptosis, MAPK signaling pathway

图片 (Image Data)



Western blot analysis of Caspase3 p12 in K562, 3T3, HeLa lysates using Cleaved-Caspase 3 p12 antibody.

注意事项 (Note)

For research use only .