

产品名称: TriMethyl-Histone H3 (Lys27) Rabbit Monoclonal Antibody
产品货号: AMRe84830



产品概述 (Summary)

产品名称 (Production Name)	TriMethyl-Histone H3 (Lys27) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TriMethyl-Histone H3 (Lys27)
别名 (Alternative Names)	H3K27me3; H3 histone; HIST1H3A; Histone cluster 1; H3a
基因 ID (Gene ID)	8350.0
蛋白 ID (SwissProt ID)	P68431.A synthetic methyl-peptide corresponding to residues surrounding Lys27 of human Histone H3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

产品名称: TriMethyl-Histone H3 (Lys27) Rabbit Monoclonal Antibody
产品货号: AMRe84830




研究背景 (Background)


H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Histone H3 (tri methyl-Lys27) in C6, Jurkat lysates using TriMethyl-Histone H3 (Lys27) antibody.



Immunocytochemistry analysis of TriMethyl-Histone H3 (Lys27) (green) in HeLa using TriMethyl-Histone H3 (Lys27) antibody, and DAPI (blue)

产品名称: TriMethyl-Histone H3 (Lys27) Rabbit Monoclonal
Antibody
产品货号: AMRe84830



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using TriMethyl-Histone H3 (Lys27) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .