

产品名称: Phospho-RSK4 (Ser232) Rabbit Monoclonal Antibody
产品货号: AMRe84835



产品概述 (Summary)

产品名称 (Production Name)	Phospho-RSK4 (Ser232) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-RSK4 (Ser232) RPS6KA6; RSK4; Ribosomal protein S6 kinase alpha-6; S6K-alpha-6; 90 kDa
别名 (Alternative Names)	ribosomal protein S6 kinase 6; p90-RSK 6; p90RSK6; Ribosomal S6 kinase 4; RSK-4; pp90RSK4
基因 ID (Gene ID)	27330.0
蛋白 ID (SwissProt ID)	Q9UK32.A synthetic phosphopeptide corresponding to residues surrounding Ser235 of human RSK4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
----------------------	---------------------------------

产品名称: Phospho-RSK4 (Ser232) Rabbit Monoclonal Antibody
产品货号: AMRe84835




蛋白分子量 (Molecular Weight) Calculated MW: 84 kDa; Observed MW: 84 kDa

研究背景 (Background)



Serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the transcription factor CREB.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-RSK4 (Ser232) in HeLa lysates using Phospho-RSK4 (Ser232) antibody.



Immunocytochemistry analysis of Phospho-RSK4 (Ser232) (green) in LNCaP using Phospho-RSK4 (Ser232) antibody, and DAPI(blue).

注意事项 (Note)

For research use only .