

产品名称: Phospho-Nrf2 (Ser40) Rabbit Monoclonal Antibody
产品货号: AMRe84837



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Nrf2 (Ser40) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Nrf2 (Ser40)
别名 (Alternative Names)	NFE2L2; NRF2; Nuclear factor erythroid 2-related factor 2; NF-E2-related factor 2; NFE2-related factor 2; HEBP1; Nuclear factor; erythroid derived 2; like 2
基因 ID (Gene ID)	4780.0
蛋白 ID (SwissProt ID)	Q16236.A synthetic phosphopeptide corresponding to residues surrounding Ser40 of human Nrf2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 100 kDa

产品名称: Phospho-Nrf2 (Ser40) Rabbit Monoclonal Antibody
产品货号: AMRe84837




研究背景 (Background)


Transcription activator that binds to antioxidant response (ARE) elements in the promoter regions of target genes. Important for the coordinated up-regulation of genes in response to oxidative stress. May be involved in the transcriptional activation of genes of the beta-globin cluster by mediating enhancer activity of hypersensitive site 2 of the beta-globin locus control region.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-Nrf2 (Ser40) in K562, HeLa lysates using Phospho-Nrf2 (Ser40) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Nrf2 (Phospho- S40) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .