

产品名称: Phospho-S6K2 (Ser371) Rabbit Monoclonal Antibody
产品货号: AMRe84850



产品概述 (Summary)

产品名称 (Production Name)	Phospho-S6K2 (Ser371) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-S6K2 (Ser371) RPS6KB2; STK14B; Ribosomal protein S6 kinase beta-2; S6K-beta-2; S6K2; 70 kDa ribosomal protein S6 kinase 2; P70S6K2; p70-S6K 2; S6 kinase-related kinase; SRK; Serine/threonine-protein kinase 14B; p70 ribosomal S6 kinase beta; S6K-beta; p70
别名 (Alternative Names)	
基因 ID (Gene ID)	6199.0
蛋白 ID (SwissProt ID)	Q9UBS0.A synthetic phosphopeptide corresponding to residues surrounding Ser371 of human P70 S6 Kinase beta

产品应用 (Application)

产品名称: Phospho-S6K2 (Ser371) Rabbit Monoclonal Antibody
产品货号: AMRe84850



稀释比 (Dilution Ratio)

WB 1:500-1:1000

蛋白分子量 (Molecular Weight)

Calculated MW: 53 kDa; Observed MW: 60 kDa


研究背景 (Background)

Phosphorylates specifically ribosomal protein S6.

研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



Western blot analysis of Phospho-p70 S6 Kinase beta/SRK (Ser371) in K562, rat Brain, C6, 3T3 lysates using Phospho-S6K2 (Ser371) antibody.

注意事项 (Note)

For research use only .