

**产品名称: Phospho-Acetyl-Coenzyme A Carboxylase (Ser79)
Rabbit Monoclonal Antibody**
产品货号: AMRe84851



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Acetyl-Coenzyme A Carboxylase (Ser79) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated,Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Acetyl Coenzyme A Carboxylase (Ser79)
别名 (Alternative Names)	ACACA; ACAC; ACC1; ACCA; Acetyl-CoA carboxylase 1; ACC1; ACC-alpha
基因 ID (Gene ID)	31.0
蛋白 ID (SwissProt ID)	Q13085.A synthetic phosphopeptide corresponding to residues surrounding Ser79 of human Acetyl Coenzyme A Carboxylase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 266 kDa; Observed MW: 277 kDa

**产品名称: Phospho-Acetyl-Coenzyme A Carboxylase (Ser79)
Rabbit Monoclonal Antibody**
产品货号: AMRe84851




研究背景 (Background)

Catalyzes the rate-limiting reaction in the biogenesis of long-chain fatty acids. Carries out three functions: biotin carboxyl carrier protein, biotin carboxylase and carboxyltransferase.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Acetyl Coenzyme A Carboxylase (Phospho-Ser79) in K562 lysates using Phospho-Acetyl Coenzyme A Carboxylase (Ser79) antibody.

注意事项 (Note)

For research use only .