

**产品名称: Phospho-Glutamate Receptor 1 (AMPA Subtype)
(Ser845) Rabbit Monoclonal Antibody**
产品货号: AMRe84852



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Glutamate Receptor 1 (AMPA Subtype) (Ser845) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Glutamate Receptor 1 (AMPA Subtype) (Ser845)
别名 (Alternative Names)	GRIA1; GLUH1; GLUR1; Glutamate receptor 1; GluR-1; AMPA-selective glutamate receptor 1; GluR-A; GluR-K1; Glutamate receptor ionotropic; AMPA 1; GluA1
基因 ID (Gene ID)	2890.0
蛋白 ID (SwissProt ID)	P42261.A synthetic phosphopeptide corresponding to residues surrounding Ser845 of human AMPA Receptor 1 (GluA1)

产品应用 (Application)

产品名称: Phospho-Glutamate Receptor 1 (AMPA Subtype)
(Ser845) Rabbit Monoclonal Antibody
产品货号: AMRe84852



稀释比 (Dilution Ratio)

WB 1:500-1:1000

蛋白分子量 (Molecular Weight)


Calculated MW: 102 kDa; Observed MW: 102 kDa

研究背景 (Background)


AMPA- (α -amino-3-hydroxy-5-methyl-4-isoxazolepropionic acid), kainate-, and NMDA- (N-methyl-D-aspartate) receptors are the three main families of ionotropic glutamate-gated ion channels. AMPA receptors (AMPA receptors) are comprised of four subunits (GluR 1-4), which assemble as homo- or hetero-tetramers to mediate the majority of fast excitatory transmissions in the central nervous system. AMPARs are implicated in synapse formation, stabilization, and plasticity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-Glutamate Receptor 1 (AMPA subtype) (Ser845) in U87-MG lysates using Phospho-Glutamate Receptor 1 (AMPA Subtype) (Ser845) antibody.



Western blot analysis of Phospho-Glutamate Receptor 1 (AMPA subtype) (Ser845) in C6 lysates using Phospho-Glutamate Receptor 1 (AMPA subtype) (Ser845) antibody.

**产品名称: Phospho-Glutamate Receptor 1 (AMPA Subtype)
(Ser845) Rabbit Monoclonal Antibody**
产品货号: AMRe84852



注意事项 (Note)

For research use only .