

产品名称: Phospho-AP2M1 (Thr156) Rabbit Monoclonal Antibody
产品货号: AMRe84853



产品概述 (Summary)

产品名称 (Production Name)	Phospho-AP2M1 (Thr156) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-AP2M1 (Thr156)
别名 (Alternative Names)	AP-2 complex subunit mu; AP-2 mu chain; Adapter-related protein complex 2 subunit mu; Adaptin-mu2; Clathrin assembly protein complex 2 mu medium chain; Clathrin coat assembly protein AP50
基因 ID (Gene ID)	1173.0
蛋白 ID (SwissProt ID)	Q96CW1.A synthetic phosphopeptide corresponding to residues surrounding Thr156 of human AP2M1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
----------------------	-----------------

产品名称: Phospho-AP2M1 (Thr156) Rabbit Monoclonal Antibody
产品货号: AMRe84853




蛋白分子量 (Molecular Weight) Calculated MW: 50 kDa; Observed MW: 50 kDa

研究背景 (Background)

Adaptins are heterotetrameric subunits of adaptors, which are complexes involved in the formation of Clathrin-coated pits for vesicle-mediated endocytosis. Involved in clathrin-dependent endocytosis in which cargo proteins are incorporated into vesicles surrounded by clathrin (clathrin-coated vesicles, CCVs) which are destined for fusion with the early endosome.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-AP2M1 (Thr156) in K562 lysates using Phospho-AP2M1 (Thr156) antibody.

注意事项 (Note)

For research use only .