

产品名称: Phospho-CD18 (Ser745) Rabbit Monoclonal Antibody
产品货号: AMRe84856



产品概述 (Summary)

产品名称 (Production Name)	Phospho-CD18 (Ser745) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-CD18 (Ser745)
别名 (Alternative Names)	ITGB2; CD18; MFI7; Integrin beta-2; Cell surface adhesion glycoproteins LFA-1/CR3/p150; 95 subunit beta; Complement receptor C3 subunit beta; CD antigen CD18
基因 ID (Gene ID)	3689.0
蛋白 ID (SwissProt ID)	P05107.A synthetic phosphopeptide corresponding to residues surrounding Ser745 of human CD18

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
----------------------	------------------------------

产品名称: Phospho-CD18 (Ser745) Rabbit Monoclonal Antibody
产品货号: AMRe84856




蛋白分子量 (Molecular Weight) Calculated MW: 85 kDa; Observed MW: 85 kDa

研究背景 (Background)

Integrin alpha-X/beta-2 recognizes the sequence G-P-R in fibrinogen alpha-chain. Integrin alpha-M/beta-2 recognizes P1 and P2 peptides of fibrinogen gamma chain. Integrin alpha-M/beta-2 is also a receptor for factor X. Integrin alpha-D/beta-2 is a receptor for ICAM3 and VCAM1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CD18 (Phospho-Ser745) in K562, rat Brain, C6, 3T3 lysates using Phospho-CD18 (Ser745) antibody.

注意事项 (Note)

For research use only .