

产品名称: Phospho-Dishevelled 2 (Ser143) Rabbit Monoclonal Antibody
产品货号: AMRe84862



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Dishevelled 2 (Ser143) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Dishevelled 2 (Ser143)
别名 (Alternative Names)	DSH homolog 2
基因 ID (Gene ID)	1856.0
蛋白 ID (SwissProt ID)	O14641.A synthetic phosphopeptide corresponding to residues surrounding Ser143 of human Dishevelled 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 79 kDa; Observed MW: 90-95 kDa

产品名称: Phospho-Dishevelled 2 (Ser143) Rabbit Monoclonal Antibody
产品货号: AMRe84862


研究背景 (Background)

Plays a role in the signal transduction pathways mediated by multiple Wnt genes. Participates both in canonical and non-canonical Wnt signaling by binding to the cytoplasmic C-terminus of frizzled family members and transducing the Wnt signal to down-stream effectors. Promotes internalization and degradation of frizzled proteins upon Wnt signaling.

研究领域 (Research Area)

Wnt signaling pathway

图片 (Image Data)



Western blot analysis of Dishevelled 2 (Phospho-Ser143) in K562 lysates using Phospho-Dishevelled 2 (Ser143) antibody.

注意事项 (Note)

For research use only .