

产品名称: Phospho-IKB alpha (Ser32) Rabbit Monoclonal Antibody
产品货号: AMRe84865



产品概述 (Summary)

产品名称 (Production Name)	Phospho-IKB alpha (Ser32) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-IKB alpha (Ser32)
别名 (Alternative Names)	NFKBIA; IKBA; MAD3; NFKBI; NF-kappa-B inhibitor alpha; I-kappa-B-alpha; Ikb-alpha; IkappaBalpha; Major histocompatibility complex enhancer-binding protein MAD3
基因 ID (Gene ID)	4792.0
蛋白 ID (SwissProt ID)	P25963.A synthetic phosphopeptide corresponding to residues surrounding Ser32 of human IKB alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
----------------------	-------------------------------

产品名称: Phospho-IKB alpha (Ser32) Rabbit Monoclonal Antibody
产品货号: AMRe84865




蛋白分子量 (Molecular Weight) Calculated MW: 36 kDa; Observed MW: 36 kDa

研究背景 (Background)

NFKB1 (MIM 164011) or NFKB2 (MIM 164012) is bound to REL (MIM 164910), RELA (MIM 164014), or RELB (MIM 604758) to form the NFKB complex. The NFKB complex is inhibited by I-kappa-B proteins (NFKBIA or NFKBIB, MIM 604495), which inactivate NF-kappa-B by trapping it in the cytoplasm.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IKB alpha (Phospho-Ser32) in 3T3, HeLa lysates using Phospho-IKB alpha (Ser32) antibody.

注意事项 (Note)

For research use only .