

产品名称: Phospho-NFAT2 (Ser237) Rabbit Monoclonal Antibody
产品货号: AMRe84869



产品概述 (Summary)

产品名称 (Production Name)	Phospho-NFAT2 (Ser237) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-NFAT2 (Ser237)
别名 (Alternative Names)	NFATC1; NFAT2; NFATC; Nuclear factor of activated T-cells; cytoplasmic 1; NF-ATc1; NFATc1; NFAT transcription complex cytosolic component; NF-ATc; NFATc
基因 ID (Gene ID)	4772.0
蛋白 ID (SwissProt ID)	O95644.A synthetic phosphopeptide corresponding to residues surrounding Ser237 of human NFAT2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
----------------------	-----------------

产品名称: Phospho-NFAT2 (Ser237) Rabbit Monoclonal Antibody
产品货号: AMRe84869



蛋白分子量 (Molecular Weight) Calculated MW: 101 kDa; Observed MW: 101 kDa


研究背景 (Background)

NFAT proteins are predominantly expressed in cells of the immune system, but are also expressed in skeletal muscle, keratinocytes, and adipocytes, regulating cell differentiation programs in these cells.

研究领域 (Research Area)

MAPK signaling pathway

图片 (Image Data)



Western blot analysis of NFAT2 (Phospho-Ser237) in K562, 3T3 lysates using Phospho-NFAT2 (Ser237) antibody.

注意事项 (Note)

For research use only .