

产品名称: Phospho-Parvin alpha (Ser8) Rabbit Monoclonal Antibody
产品货号: AMRe84872



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Parvin alpha (Ser8) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Parvin alpha (Ser8)
别名 (Alternative Names)	MXRA2; CH-ILKBP
基因 ID (Gene ID)	55742.0
蛋白 ID (SwissProt ID)	Q9NVD7.A synthetic phosphopeptide corresponding to residues surrounding Ser8 of human Parvin alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 42 kDa; Observed MW: 42 kDa


产品名称: Phospho-Parvin alpha (Ser8) Rabbit Monoclonal Antibody
产品货号: AMRe84872

研究背景 (Background)

Plays a role in sarcomere organization and in smooth muscle cell contraction. Required for normal development of the embryonic cardiovascular system, and for normal septation of the heart outflow tract. Plays a role in sprouting angiogenesis and is required for normal adhesion of vascular smooth muscle cells to endothelial cells during blood vessel development . Plays a role in the reorganization of the actin cytoskeleton, formation of lamellipodia and ciliogenesis. Plays a role in the establishment of cell polarity, cell adhesion, cell spreading, and directed cell migration.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Parvin alpha (Phospho-Ser8) in U2OS, mouse brain, rat brain, PC-12, A549 lysates using Phospho-Parvin alpha (Ser8) antibody.

注意事项 (Note)

For research use only .