

产品名称: Phospho-PKC gamma (Thr514) Rabbit Monoclonal Antibody
产品货号: AMRe84874



产品概述 (Summary)

产品名称 (Production Name)	Phospho-PKC gamma (Thr514) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-PKC gamma (Thr514)
别名 (Alternative Names)	PRKCG; PKCG; Protein kinase C gamma type; PKC-gamma
基因 ID (Gene ID)	5582.0
蛋白 ID (SwissProt ID)	P05129.A synthetic phosphopeptide corresponding to residues surrounding Thr514 of human PKC

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 78 kDa; Observed MW: 78-85 kDa

产品名称: Phospho-PKC gamma (Thr514) Rabbit Monoclonal Antibody
产品货号: AMRe84874




研究背景 (Background)

This is a calcium-activated, phospholipid-dependent, serine- and threonine-specific enzyme.


研究领域 (Research Area)

MAPK signaling pathway


图片 (Image Data)



Western blot analysis of Phospho-PKC gamma (Thr514) in mouse brain lysates using Phospho-PKC gamma (Thr514) antibody.



Western blot analysis of Phospho-PKC gamma (Thr514) in C6, HeLa lysates using Phospho-PKC gamma (Thr514) antibody.



产品名称: Phospho-PKC gamma (Thr514) Rabbit Monoclonal Antibody
产品货号: AMRe84874



Immunohistochemistry analysis of paraffin-embedded Human Brain using PKC (Phospho-Thr514) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .