

产品名称: Phospho-Progesterone Receptor (Ser190) Rabbit Monoclonal Antibody
产品货号: AMRe84878



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Progesterone Receptor (Ser190) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Progesterone Receptor (Ser190)
别名 (Alternative Names)	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3
基因 ID (Gene ID)	5241.0
蛋白 ID (SwissProt ID)	P06401.A synthetic phosphopeptide corresponding to residues surrounding Ser190 of human Progesterone Receptor

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 99 kDa; Observed MW: 90, 118 kDa

产品名称: Phospho-Progesterone Receptor (Ser190) Rabbit
Monoclonal Antibody
产品货号: AMRe84878




研究背景 (Background)

Progesterone plays a central role in the reproductive events associated with the establishment and maintenance of pregnancy. Progesterone receptor, a member of the steroid receptor superfamily, mediates the physiologic effects of progesterone.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Progesterone Receptor (Phospho-Ser190) in T47D lysates using Phospho-Progesterone Receptor (Ser190) antibody.

注意事项 (Note)

For research use only .