

产品名称: Phospho-Vimentin (Ser56) Rabbit Monoclonal Antibody
产品货号: AMRe84884



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Vimentin (Ser56) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Vimentin (Ser56)
别名 (Alternative Names)	VIM; Vimentin
基因 ID (Gene ID)	7431.0
蛋白 ID (SwissProt ID)	P08670.A synthetic phosphopeptide corresponding to residues surrounding Ser56 of human Vimentin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 54 kDa

产品名称: Phospho-Vimentin (Ser56) Rabbit Monoclonal Antibody
产品货号: AMRe84884




研究背景 (Background)


Vimentin an intermediate filament protein. Intermediate filament proteins are expressed in a tissue-specific manner. Desmin is the subunit specific for muscle and vimentin the subunit specific for mesenchymal tissue.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-Vimentin (Ser56) in U251MG lysates using Phospho-Vimentin (Ser56) antibody.



Immunocytochemistry analysis of Phospho-Vimentin (Ser56) (green) in Hela using Phospho-Vimentin (Ser56) antibody, and DAPI(blue).

注意事项 (Note)

For research use only .