

产品名称: Phospho-Cyclin B1 (Ser126) Rabbit Monoclonal Antibody
产品货号: AMRe84894



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Cyclin B1 (Ser126) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Cyclin B1 (Ser126)
别名 (Alternative Names)	CCNB1; CCNB; G2/mitotic-specific cyclin-B1
基因 ID (Gene ID)	891.0
蛋白 ID (SwissProt ID)	P14635.A synthetic phosphopeptide corresponding to residues surrounding Ser126 of human Cyclin B1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 55 kDa

产品名称: Phospho-Cyclin B1 (Ser126) Rabbit Monoclonal Antibody
产品货号: AMRe84894




研究背景 (Background)

Essential for the control of the cell cycle at the G2/M (mitosis) transition. a member of the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Cyclin B1 (Phospho-Ser126) in HeLa, Ramos lysates using Phospho-Cyclin B1 (Ser126) antibody.

注意事项 (Note)

For research use only .