

产品名称: Phospho-ErbB 4 (Tyr1162) Rabbit Monoclonal Antibody
产品货号: AMRe84899



产品概述 (Summary)

产品名称 (Production Name)	Phospho-ErbB 4 (Tyr1162) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-ErbB 4 (Tyr1162)
别名 (Alternative Names)	ERBB4; HER4; Receptor tyrosine-protein kinase erbB-4; Proto-oncogene-like protein c-ErbB-4; Tyrosine kinase-type cell surface receptor HER4; p180erbB4
基因 ID (Gene ID)	2066.0
蛋白 ID (SwissProt ID)	Q15303.A synthetic phosphopeptide corresponding to residues surrounding Tyr1162 of human ErbB 4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 147 kDa; Observed MW: 180 kDa

产品名称: Phospho-ErbB 4 (Tyr1162) Rabbit Monoclonal Antibody
产品货号: AMRe84899




研究背景 (Background)

The HER4/ERBB4 gene is a member of the type I receptor tyrosine kinase subfamily that includes EGFR (MIM 131550), ERBB2 (MIM 164870), and ERBB3 (MIM 190151). It encodes a receptor for NDF/heregulin (MIM 142445).


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of Phospho-ErbB 4 (Tyr1162) in HeLa, A549 lysates using Phospho-ErbB 4 (Tyr1162) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using ErbB 4 (Phospho-Tyr1162) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .