

**产品名称: Phospho-Ezrin/Radixin/Moesin  
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody**  
**产品货号: AMRe84900**



## 产品概述 (Summary)

|                        |                                                                                |
|------------------------|--------------------------------------------------------------------------------|
| 产品名称 (Production Name) | Phospho-Ezrin/Radixin/Moesin (Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody                                         |
| 宿主 (Host)              | Rabbit                                                                         |
| 应用 (Application)       | WB,IHC,IP                                                                      |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat,Hamster                                                        |

## 产品性能 (Performance)

|                     |                                                                                            |
|---------------------|--------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                               |
| 修饰 (Modification)   | Phosphorylated                                                                             |
| 同种型 (Isotype)       | IgG                                                                                        |
| 克隆 (Clonality)      | Monoclonal                                                                                 |
| 形式 (Form)           | Liquid                                                                                     |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.   |
| 储存溶液 (Buffer)       | Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                      |

## 免疫原信息 (Immunogen)

|                        |                                                                                               |
|------------------------|-----------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | Phospho-Ezrin/Radixin/Moesin (Thr567/Thr564/Thr558)                                           |
| 别名 (Alternative Names) | EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81                                                   |
| 基因 ID (Gene ID)        | 7430.0                                                                                        |
| 蛋白 ID (SwissProt ID)   | P15311.A synthetic phosphopeptide corresponding to residues surrounding Thr567 of human Ezrin |

## 产品应用 (Application)

|                          |                                               |
|--------------------------|-----------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20   |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 69 kDa; Observed MW: 75,80 kDa |


产品名称: Phospho-Ezrin/Radixin/Moesin  
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody  
产品货号: AMRe84900

## 研究背景 (Background)


Probably involved in connections of major cytoskeletal structures to the plasma membrane. In epithelial cells, required for the formation of microvilli and membrane ruffles on the apical pole. Along with PLEKHG6, required for normal macropinocytosis.

## 研究领域 (Research Area)

## 图片 (Image Data)




Western blot analysis of Phospho-Ezrin (Thr567)/ Radixin (Thr564)/ Moesin (Thr558) in 3T3, 3T3+PDGF lysates using Phospho-Ezrin/Radixin/Moesin (Thr567/Thr564/Thr558) antibody.



Western blot analysis of Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) in HeLa, CHO-K1, C6 lysates using Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) antibody.

**产品名称: Phospho-Ezrin/Radixin/Moesin  
(Thr567/Thr564/Thr558) Rabbit Monoclonal Antibody**  
**产品货号: AMRe84900**

---



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Ezrin (Phospho-Thr567)/ Radixin (Phospho-Thr564)/ Moesin (Phospho-Thr558) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

#### **注意事项 (Note)**

For research use only .