

产品名称: Phospho-Filamin A/B (Ser2152/Ser2107) Rabbit Monoclonal Antibody
产品货号: AMRe84902



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Filamin A/B (Ser2152/Ser2107) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Filamin A/B (Ser2152/Ser2107)
别名 (Alternative Names)	FLN; FMD; MNS; OPD; ABPX; CSBS; CVD1; FGS2; FLN1; NHBP; OPD1; OPD2; XLVD; XMVD; FLN-A; ABP-280
基因 ID (Gene ID)	2317.0
蛋白 ID (SwissProt ID)	O75369.A synthetic phosphopeptide corresponding to residues surrounding Ser2152 of human Filamin A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 278 kDa; Observed MW: 281 kDa

产品名称: Phospho-Filamin A/B (Ser2152/Ser2107) Rabbit
Monoclonal Antibody
产品货号: AMRe84902




研究背景 (Background)

Connects cell membrane constituents to the actin cytoskeleton.


研究领域 (Research Area)

MAPK signaling pathway

图片 (Image Data)




Western blot analysis of Filamin A (Phospho-Ser2152)/Filamin B (Phospho-Ser2107) in HeLa, A549, U2OS, C6 lysates using Filamin A (Phospho-Ser2152)/Filamin B (Phospho-Ser2107) antibody.



Immunocytochemistry analysis of Phospho-Filamin A/B (Ser2152/Ser2107) (green) in A549 using Phospho-Filamin A/B (Ser2152/Ser2107) antibody, and DAPI (blue).

**产品名称: Phospho-Filamin A/B (Ser2152/Ser2107) Rabbit
Monoclonal Antibody**
产品货号: AMRe84902



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Filamin A (Phospho-Ser2152)/Filamin B (Phospho-Ser2107) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .