

产品名称: Phospho-GATA3 (Ser308) Rabbit Monoclonal Antibody
产品货号: AMRe84903



产品概述 (Summary)

产品名称 (Production Name)	Phospho-GATA3 (Ser308) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-GATA3 (Ser308)
别名 (Alternative Names)	GATA3; Trans-acting T-cell-specific transcription factor GATA-3; GATA-binding factor 3
基因 ID (Gene ID)	2625.0
蛋白 ID (SwissProt ID)	P23771.A synthetic phosphopeptide corresponding to residues surrounding Ser308 of human GATA3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 48 kDa

产品名称: Phospho-GATA3 (Ser308) Rabbit Monoclonal Antibody
产品货号: AMRe84903




研究背景 (Background)


Transcriptional activator which binds to the enhancer of the T-cell receptor alpha and delta genes. Binds to the consensus sequence 5'-AGATAG-3'. Required for the T-helper 2 (Th2) differentiation process following immune and inflammatory responses.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-GATA3 (Ser308) in A549, MCF-7, U251 lysates using Phospho-GATA3 (Ser308) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using GATA3 (Phospho-Ser308) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .