

**产品名称: Phospho-Glucocorticoid Receptor (Ser226) Rabbit Monoclonal Antibody**  
**产品货号: AMRe84904**



## 产品概述 (Summary)

|                        |                                                                     |
|------------------------|---------------------------------------------------------------------|
| 产品名称 (Production Name) | Phospho-Glucocorticoid Receptor (Ser226) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody                              |
| 宿主 (Host)              | Rabbit                                                              |
| 应用 (Application)       | WB,IHC                                                              |
| 种属反应性 (Reactivity)     | Human,Rat                                                           |

## 产品性能 (Performance)

|                     |                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                 |
| 修饰 (Modification)   | Phosphorylated                                                                               |
| 同种型 (Isotype)       | IgG                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                   |
| 形式 (Form)           | Liquid                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.     |
| 储存溶液 (Buffer)       | Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                 |
|------------------------|-----------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | Phospho-Glucocorticoid Receptor (Ser226)                                                                        |
| 别名 (Alternative Names) | NR3C1; GRL; Glucocorticoid receptor; GR; Nuclear receptor subfamily 3 group C member 1                          |
| 基因 ID (Gene ID)        | 2908.0                                                                                                          |
| 蛋白 ID (SwissProt ID)   | P04150.A synthetic phosphopeptide corresponding to residues surrounding Ser226 of human Glucocorticoid Receptor |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000, IHC 1:50-1:100            |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 86 kDa; Observed MW: 94 kDa |

产品名称: Phospho-Glucocorticoid Receptor (Ser226) Rabbit  
Monoclonal Antibody  
产品货号: AMRe84904




## 研究背景 (Background)

The protein encoded by this gene is a receptor for glucocorticoids and can act as both a transcription factor and a regulator of other transcription factors. The encoded protein can bind DNA as a homodimer or as a heterodimer with another protein such as the retinoid X receptor. This protein can also be found in heteromeric cytoplasmic complexes along with heat shock factors and immunophilins.


## 研究领域 (Research Area)

Apoptosis

## 图片 (Image Data)



Western blot analysis of Phospho-Glucocorticoid Receptor (Ser226) in HeLa, A549, C6 lysates using Phospho-Glucocorticoid Receptor (Ser226) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Glucocorticoid Receptor (Phospho-Ser226) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## 注意事项 (Note)

For research use only .

**产品名称: Phospho-Glucocorticoid Receptor (Ser226) Rabbit  
Monoclonal Antibody**  
**产品货号: AMRe84904**

---

