

产品名称: Phospho-Jun/JunD (Ser73/Ser100) Rabbit Monoclonal Antibody
产品货号: AMRe84907



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Jun/JunD (Ser73/Ser100) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Jun/JunD (Ser73/Ser100)
别名 (Alternative Names)	JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39; JUND; Transcription factor jun-D
基因 ID (Gene ID)	3725/3727
蛋白 ID (SwissProt ID)	P05412/P17535.A synthetic phosphopeptide corresponding to residues surrounding Ser73 of human Jun

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
----------------------	------------------------------

产品名称: Phospho-Jun/JunD (Ser73/Ser100) Rabbit Monoclonal Antibody
产品货号: AMRe84907




蛋白分子量 (Molecular Weight) Calculated MW: 36 kDa; Observed MW: 38,42,48 kDa

研究背景 (Background)

c-jun is a transcription factor that recognizes and binds to the enhancer heptamer motif 5'-TGA[CG]TCA-3'. Promotes activity of NR5A1 when phosphorylated by HIPK3 leading to increased steroidogenic gene expression upon cAMP signaling pathway stimulation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Jun (Phospho-Ser73)/JunD (Phospho-Ser100) in A549, HL-60, U251, U87-MG lysates using Phospho-Jun/JunD (Ser73/Ser100) antibody.

注意事项 (Note)

For research use only .