

产品名称: Phospho-c-Met (Tyr1349) Rabbit Monoclonal Antibody
产品货号: AMRe84909



产品概述 (Summary)

产品名称 (Production Name)	Phospho-c-Met (Tyr1349) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-c-Met (Tyr1349)
别名 (Alternative Names)	MET; Hepatocyte growth factor receptor; HGF receptor; HGF/SF receptor; Proto-oncogene c-Met; Scatter factor receptor; SF receptor; Tyrosine-protein kinase Met
基因 ID (Gene ID)	4233.0
蛋白 ID (SwissProt ID)	P08581.A synthetic phosphopeptide corresponding to residues surrounding Tyr1349 of human Met

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
----------------------	-------------------------------

产品名称: Phospho-c-Met (Tyr1349) Rabbit Monoclonal Antibody
产品货号: AMRe84909



蛋白分子量 (Molecular Weight) Calculated MW: 156 kDa; Observed MW: 170,140 kDa


研究背景 (Background)

The proto-oncogene MET product is the hepatocyte growth factor receptor and encodes tyrosine-kinase activity. The primary single chain precursor protein is post-translationally cleaved to produce the alpha and beta subunits, which are disulfide linked to form the mature receptor.

研究领域 (Research Area)

TGF-beta signaling pathway, PI3K-Akt signaling pathway, Hippo signaling pathway

图片 (Image Data)



Western blot analysis of Phospho-Met (c-Met) (Tyr1349) in HeLa, A549 lysates using Phospho-c-Met (Tyr1349) antibody.

注意事项 (Note)

For research use only .