

产品名称: Phospho-NDRG1 (Ser330) Rabbit Monoclonal Antibody
产品货号: AMRe84911



产品概述 (Summary)

产品名称 (Production Name)	Phospho-NDRG1 (Ser330) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-NDRG1 (Ser330)
别名 (Alternative Names)	Protein NDRG1; Differentiation-related gene 1 protein; DRG-1; RTP; Rit42; NDRG1; CAP43; DRG1; RTP; targ1; TDD5; tdds
基因 ID (Gene ID)	10397.0
蛋白 ID (SwissProt ID)	Q92597.A synthetic phosphopeptide corresponding to residues surrounding Ser330 of human NDRG1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 43 kDa

产品名称: Phospho-NDRG1 (Ser330) Rabbit Monoclonal Antibody
产品货号: AMRe84911




研究背景 (Background)


May have a growth inhibitory role.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Phospho-NDRG1 (Ser330) in mouse kidney lysates using Phospho-NDRG1 (Ser330) antibody.




Western blot analysis of Phospho-NDRG1(Ser330) in HeLa, A549, MCF-7, 239T, C2C12 lysates using Phospho-NDRG1(Ser330) antibody

产品名称: Phospho-NDRG1 (Ser330) Rabbit Monoclonal Antibody
产品货号: AMRe84911



Immunocytochemistry analysis of Phospho-NDRG1 (Ser330) (green) in 293T using Phospho-NDRG1 (Ser330) antibody, and DAPI(blue).



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using NDRG1 (Phospho-Ser330) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .