

产品名称: Phospho-Histone H3 (Ser10/Thr11) Rabbit Monoclonal Antibody
产品货号: AMRe84921



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Histone H3 (Ser10/Thr11) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Histone H3 (Ser10/Thr11)
别名 (Alternative Names)	H3 histone; family 3A; H3 histone; family 3B (H3.3B); H3.3A; H3.3B; H33; H3F3; H3F3A; H3F3B; Histone H3.3
基因 ID (Gene ID)	8350.0
蛋白 ID (SwissProt ID)	P68431.A synthetic phosphopeptide corresponding to residues surrounding Ser10/Thr11 of human Histone H3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

产品名称: Phospho-Histone H3 (Ser10/Thr11) Rabbit
Monoclonal Antibody
产品货号: AMRe84921




研究背景 (Background)


H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Phospho-Histone H3 (Ser10/Thr11) in 3T3 lysates using Phospho-Histone H3 (Ser10/Thr11)antibody.




Western blot analysis of Histone H3 (Phospho-Ser10/Thr11) in HeLa, A549, HL-60, C6 lysates using Histone H3 (Phospho-Ser10/Thr11) antibody

**产品名称: Phospho-Histone H3 (Ser10/Thr11) Rabbit
Monoclonal Antibody**
产品货号: AMRe84921



Immunocytochemistry analysis of Phospho-Histone H3 (Ser10/Thr11) (green) in HeLa using Phospho-Histone H3 (Ser10/Thr11) antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Histone H3 (Phospho-Ser10/Thr11) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .