

**产品名称: Phospho-ERK1/2
(Thr202/Tyr204)/(Thr185/Tyr187) Rabbit Monoclonal
Antibody**
产品货号: AMRe84925



产品概述 (Summary)

产品名称 (Production Name)	Phospho-ERK1/2 (Thr202/Tyr204)/(Thr185/Tyr187) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-ERK1/2 (Thr202/Tyr204)/(Thr185/Tyr187)
别名 (Alternative Names)	MAPK1/MAPK3
基因 ID (Gene ID)	5595/5594
蛋白 ID (SwissProt ID)	P27361/P28482.A synthetic phosphopeptide corresponding to residues surrounding Thr202/Tyr204 of human Erk1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
----------------------	------------------------------

**产品名称: Phospho-ERK1/2
(Thr202/Tyr204)/(Thr185/Tyr187) Rabbit Monoclonal
Antibody**
产品货号: AMRe84925


蛋白分子量 (Molecular Weight) Calculated MW: 44,42 kDa; Observed MW: 44,42 kDa

研究背景 (Background)

Serine/threonine kinase which acts as an essential component of the MAP kinase signal transduction pathway. MAPK1/ERK2 and MAPK3/ERK1 are the 2 MAPKs which play an important role in the MAPK/ERK cascade. They participate also in a signaling cascade initiated by activated KIT and KITLG/SCF. Depending on the cellular context, the MAPK/ERK cascade mediates diverse biological functions such as cell growth, adhesion, survival and differentiation through the regulation of transcription, translation, cytoskeletal rearrangements.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Erk1 (pT202/pY204)/Erk2 (pT185/pY187) in C6, HeLa lysates using Phospho-ERK1/2 (Thr202/Tyr204)/(Thr185/Tyr187) antibody.

注意事项 (Note)

For research use only .