

产品名称: Phospho-Histone H2A.X (Ser139) Rabbit Monoclonal Antibody
产品货号: AMRe84926



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Histone H2A.X (Ser139) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Histone H2A.X (Ser139)
别名 (Alternative Names)	H2A.X; H2AFX; H2a/x; HIST5-2AX; Histone H2A.X; gamma H2A.X
基因 ID (Gene ID)	3014.0
蛋白 ID (SwissProt ID)	P16104.A synthetic phosphopeptide corresponding to residues surrounding Ser139 of human gamma H2A.X

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

产品名称: Phospho-Histone H2A.X (Ser139) Rabbit Monoclonal Antibody
产品货号: AMRe84926




研究背景 (Background)


Variant histone H2A which replaces conventional H2A in a subset of nucleosomes. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of gamma H2A.X (Phospho-Ser139) in HepG2 lysates using Phospho-Histone H2A.X (Ser139) antibody.



Western blot analysis of gamma H2A.X (Phospho-Ser139) in Jurkat, THP-1, C6, 3T3 lysates using gamma H2A.X (Phospho-Ser139) antibody.

注意事项 (Note)

For research use only .