

产品名称: Phospho-TAK1 (Ser439) Rabbit Monoclonal Antibody
产品货号: AMRe84928



产品概述 (Summary)

产品名称 (Production Name)	Phospho-TAK1 (Ser439) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-TAK1 (Ser439)
别名 (Alternative Names)	MAP3K7; TAK1; Mitogen-activated protein kinase kinase kinase 7; Transforming growth factor-beta-activated kinase 1; TGF-beta-activated kinase 1
基因 ID (Gene ID)	6885.0
蛋白 ID (SwissProt ID)	O43318.A synthetic phosphopeptide corresponding to residues surrounding Ser439 of human TAK1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
----------------------	---

产品名称: Phospho-TAK1 (Ser439) Rabbit Monoclonal Antibody
产品货号: AMRe84928



蛋白分子量 (Molecular Weight) Calculated MW: 67 kDa; Observed MW: 78 kDa


研究背景 (Background)

Component of a protein kinase signal transduction cascade. Mediator of TRAF6 and TGF-beta signal transduction. Activates IKKB and MAPK8 in response to TRAF6 signaling. Stimulates NF-kappa-B activation and the p38 MAPK pathway. In osmotic stress signaling, plays a major role in the activation of MAPK8/JNK, but not that of NF-kappa-B.

研究领域 (Research Area)

Apoptosis, TGF-beta signaling pathway, MAPK signaling pathway

图片 (Image Data)



Western blot analysis of TAK1 (Phospho-Ser439) in 3T3, C6, MCF-7 lysates using Phospho-TAK1 (Ser439) antibody.

注意事项 (Note)

For research use only .