

产品名称: Phospho-PKC (Ser660) Rabbit Monoclonal Antibody
产品货号: AMRe84931



产品概述 (Summary)

产品名称 (Production Name)	Phospho-PKC (Ser660) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-PKC (Ser660)
别名 (Alternative Names)	PRKCA; PKCA; PRKACA; Protein kinase C alpha type; PKC-A; PKC-alpha
基因 ID (Gene ID)	5578/5579/5580/5581/5583/5588
蛋白 ID (SwissProt ID)	P17252/P05771/Q05655/Q02156/P24723/Q04759.A synthetic phosphopeptide corresponding to residues surrounding Ser660 of human PKC

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa; Observed MW: 78-85 kDa

产品名称: Phospho-PKC (Ser660) Rabbit Monoclonal Antibody
产品货号: AMRe84931




研究背景 (Background)


PKC alpha is an AGC kinase of the PKC family. A classical PKC downstream of many mitogenic and receptors. Classical PKCs are calcium-dependent enzymes that are activated by phosphatidylserine, diacylglycerol and phorbol esters.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PKC (Phospho-Ser660) in HeLa, A549, HL-60, U2OS, C6 lysates using Phospho-PKC (Ser660) antibody.



Immunohistochemistry analysis of paraffin-embedded Human Brain using PKC (Phospho-Ser660) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .