

产品名称: UFM1 Rabbit Monoclonal Antibody
产品货号: AMRe85109



产品概述 (Summary)

产品名称 (Production Name)	UFM1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UFM1
别名 (Alternative Names)	Ubiquitin-fold modifier 1; C13orf20
基因 ID (Gene ID)	51569.0
蛋白 ID (SwissProt ID)	P61960.A synthetic peptide of human UFM1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 9 kDa; Observed MW: 9 kDa

研究背景 (Background)

Ubiquitin-like modifier which can be covalently attached via an isopeptide bond to lysine residues of substrate proteins as


产品名称: UFM1 Rabbit Monoclonal Antibody
产品货号: AMRe85109




a monomer or a lysine-linked polymer.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of UFM1 in 3T3, Raw264.7 lysates using UFM1 antibody.



Western blot analysis of UFM1 in Jurkat lysates using UFM1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using UFM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: UFM1 Rabbit Monoclonal Antibody
产品货号: AMRe85109



注意事项 (Note)

For research use only .