

产品名称: ZNF364 Rabbit Monoclonal Antibody
产品货号: AMRe85110



产品概述 (Summary)

产品名称 (Production Name)	ZNF364 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF364
别名 (Alternative Names)	E3 ubiquitin-protein ligase RNF115; Rabring 7; Zinc finger protein 364
基因 ID (Gene ID)	27246.0
蛋白 ID (SwissProt ID)	Q9Y4L5.A synthetic peptide of human ZNF364

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 34 kDa

研究背景 (Background)

E3 ubiquitin-protein ligase that mediates E2-dependent, 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of substrates


产品名称: ZNF364 Rabbit Monoclonal Antibody
产品货号: AMRe85110




and may play a role in diverse biological processes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ZNF364/RFN115 in PC-12 lysates using ZNF364 antibody.



Western blot analysis of ZNF364/RFN115 in Jurkat, Hela lysates using ZNF364/RFN115 antibody.

注意事项 (Note)

For research use only .