

产品名称: NAT10 Rabbit Monoclonal Antibody
产品货号: AMRe85116



产品概述 (Summary)

产品名称 (Production Name)	NAT10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NAT10
别名 (Alternative Names)	ALP; hALP; N acetyltransferase like protein; NAT10; NET43
基因 ID (Gene ID)	55226.0
蛋白 ID (SwissProt ID)	Q9H0A0.Recombinant protein of human NAT10

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 116 kDa; Observed MW: 116 kDa

研究背景 (Background)

Has protein acetyltransferase activity in vitro. Can acetylate both histones and microtubules. Histone acetylation may


产品名称: NAT10 Rabbit Monoclonal Antibody
产品货号: AMRe85116




regulate transcription and mitotic chromosome de-condensation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NAT10 in Jurkat, C6, HeLa lysates using NAT10 antibody.



Immunocytochemistry analysis of NAT10(green) in HeLa using NAT10 antibody, and DAPI(blue)

注意事项 (Note)

For research use only .