

产品名称: TNFAIP8 Rabbit Monoclonal Antibody
产品货号: AMRe85118



产品概述 (Summary)

产品名称 (Production Name)	TNFAIP8 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNFAIP8
别名 (Alternative Names)	TNF alpha-induced protein 8; NDED; SCC-S2
基因 ID (Gene ID)	25816.0
蛋白 ID (SwissProt ID)	O95379.A synthetic peptide of human TNFAIP8

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 21-23 kDa

研究背景 (Background)

Acts as a negative mediator of apoptosis and may play a role in tumor progression.


产品名称: TNFAIP8 Rabbit Monoclonal Antibody
产品货号: AMRe85118



研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of TNFAIP8 in CEM, K562 lysates using TNFAIP8 antibody.

注意事项 (Note)

For research use only .