

产品名称: Importin 9 Rabbit Monoclonal Antibody
产品货号: AMRe85125



产品概述 (Summary)

产品名称 (Production Name)	Importin 9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Importin 9
别名 (Alternative Names)	importin 9; Imp9; Ran-binding protein 9; RanBP9
基因 ID (Gene ID)	55705.0
蛋白 ID (SwissProt ID)	Q96P70.A synthetic peptide of human Importin 9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 116 kDa; Observed MW: 116 kDa

研究背景 (Background)


Nuclear transport receptor that mediates nuclear import of proteins, such as histones, proteasome and actin.

产品名称: Importin 9 Rabbit Monoclonal Antibody
产品货号: AMRe85125




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Importin 9/RANBP9 in Jurkat, HeLa lysates using Importin 9/RANBP9 antibody.



Immunocytochemistry analysis of Importin 9 (green) in A549 using Importin 9 antibody, and DAPI (blue).

注意事项 (Note)

For research use only.