

产品名称: MCM5 Rabbit Monoclonal Antibody
产品货号: AMRe85126



产品概述 (Summary)

产品名称 (Production Name)	MCM5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MCM5
别名 (Alternative Names)	CDC46; MCM5; MGC5315; P1 CDC46
基因 ID (Gene ID)	4174.0
蛋白 ID (SwissProt ID)	P33992.A synthetic peptide of human MCM5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 82 kDa; Observed MW: 82 kDa

研究背景 (Background)


MCM5 Belongs to the MCM family.

产品名称: MCM5 Rabbit Monoclonal Antibody
产品货号: AMRe85126




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of MCM5 in Jurkat, C6, HeLa lysates using MCM5 antibody.



Immunocytochemistry analysis of MCM5(green) in HeLa using MCM5 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: MCM5 Rabbit Monoclonal Antibody
产品货号: AMRe85126



注意事项 (Note)

For research use only .