

产品名称: MSP Rabbit Monoclonal Antibody
产品货号: AMRe85133



产品概述 (Summary)

产品名称 (Production Name)	MSP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MSP
别名 (Alternative Names)	serine threonine kinase 3/4; KRS1; MST2/KRS2; MST1; YSK3; TIAC
基因 ID (Gene ID)	6789.0
蛋白 ID (SwissProt ID)	Q13043.A synthetic peptide of human Serine/threonine-protein kinase 4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa; Observed MW: 60 kDa

研究背景 (Background)

Stress-activated, pro-apoptotic kinase which, following caspase-cleavage, enters the nucleus and induces chromatin

产品名称: MSP Rabbit Monoclonal Antibody
产品货号: AMRe85133




condensation followed by internucleosomal DNA fragmentation.


研究领域 (Research Area)

Apoptosis, MAPK signaling pathway

图片 (Image Data)



Western blot analysis of MSP/MST1 in K562, rat Brain, C6, 3T3, HeLa lysates using MSP/MST1 antibody.



Immunocytochemistry analysis of MSP (green) in K562 using MSP antibody, and DAPI (blue).

注意事项 (Note)

For research use only .