

产品名称: ATG13 Rabbit Monoclonal Antibody
产品货号: AMRe85134



产品概述 (Summary)

产品名称 (Production Name)	ATG13 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATG13
别名 (Alternative Names)	Autophagy-related protein 13
基因 ID (Gene ID)	9776.0
蛋白 ID (SwissProt ID)	O75143.Recombinant protein of human KIAA0652

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa; Observed MW: 72 kDa

研究背景 (Background)

Autophagy factor required for autophagosome formation and mitophagy.


产品名称: ATG13 Rabbit Monoclonal Antibody
产品货号: AMRe85134




研究领域 (Research Area)

Autophagy, mTOR signaling pathway

图片 (Image Data)



Western blot analysis of ATG13 in Raw264.7 lysates using ATG13 antibody.



Western blot analysis of ATG13 in K562, C6, HeLa lysates using ATG13 antibody.

注意事项 (Note)

For research use only.