

产品名称: Asparagine Synthetase Rabbit Monoclonal Antibody
产品货号: AMRe85138



产品概述 (Summary)

产品名称 (Production Name)	Asparagine Synthetase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Asparagine Synthetase
别名 (Alternative Names)	TS11; ASNSD; ASNS; Asparagine synthetase [glutamine-hydrolyzing]
基因 ID (Gene ID)	440.0
蛋白 ID (SwissProt ID)	P08243.A synthetic peptide of human Asparagine synthetase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 64 kDa; Observed MW: 64 kDa

研究背景 (Background)


产品名称: Asparagine Synthetase Rabbit Monoclonal Antibody
产品货号: AMRe85138




The protein encoded by this gene is involved in the synthesis of asparagine. This gene complements a mutation in the temperature-sensitive hamster mutant ts11, which blocks progression through the G1 phase of the cell cycle at nonpermissive temperature. Alternatively spliced transcript variants have been described for this gene.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Asparagine Synthetase in rat pancreas lysates using Asparagine Synthetase antibody.



Western blot analysis of Asparagine Synthetase in K562, HeLa lysates using Asparagine Synthetase antibody.



产品名称: Asparagine Synthetase Rabbit Monoclonal Antibody
产品货号: AMRe85138



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Asparagine synthetase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .