

产品名称: Menin Rabbit Monoclonal Antibody
产品货号: AMRe85144



产品概述 (Summary)

产品名称 (Production Name)	Menin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Menin
别名 (Alternative Names)	MEN1; SCG2; Menin
基因 ID (Gene ID)	4221.0
蛋白 ID (SwissProt ID)	O00255.A synthetic peptide of human Menin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 68 kDa

研究背景 (Background)

Menin, the protein product of the MEN1 gene, is a component of the mixed-lineage leukemia protein (MLL)-containing


产品名称: Menin Rabbit Monoclonal Antibody
产品货号: AMRe85144



histone methyltransferase complex that facilitates methylation of histone H3 Lys4 to promote transcriptional activation. Menin functions to suppress proliferation of pancreatic islet cells, at least in part through MLL-mediated activation of the p18INK4c (p18) and p27CIP/KIP (p27) cyclin-dependent kinase inhibitor genes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Menin in K562, C6 lysates using Menin antibody.

注意事项 (Note)

For research use only .