

产品名称: P2Y6 Rabbit Monoclonal Antibody
产品货号: AMRe85149



产品概述 (Summary)

产品名称 (Production Name)	P2Y6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	P2Y6
别名 (Alternative Names)	P2ry6; P2Y6; P2Y6 receptor; PP2891
基因 ID (Gene ID)	5031.0
蛋白 ID (SwissProt ID)	Q15077.A synthetic peptide of human P2Y6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 36 kDa; Observed MW: 36 kDa

研究背景 (Background)

Receptor for extracellular UDP > UTP > ATP. The activity of this receptor is mediated by G proteins which activate a


产品名称: P2Y6 Rabbit Monoclonal Antibody
产品货号: AMRe85149




phosphatidylinositol-calcium second messenger system.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of P2Y6 in K562, HeLa lysates using P2Y6 antibody.



Immunocytochemistry analysis of P2Y6 (green) in A549 using P2Y6 antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using P2Y6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: P2Y6 Rabbit Monoclonal Antibody
产品货号: AMRe85149



注意事项 (Note)

For research use only .