

产品名称: SNX5 Rabbit Monoclonal Antibody
产品货号: AMRe85158



产品概述 (Summary)

产品名称 (Production Name)	SNX5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNX5
别名 (Alternative Names)	Sorting nexin-5
基因 ID (Gene ID)	27131.0
蛋白 ID (SwissProt ID)	Q9Y5X3.Recombinant protein of human SNX5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)


Involved in several stages of intracellular trafficking.

产品名称: SNX5 Rabbit Monoclonal Antibody
产品货号: AMRe85158




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of SNX5 in mouse spleen, mouse brain lysates using SNX5 antibody.



Western blot analysis of SNX5 in K562, rat Brain, C6, HeLa lysates using SNX5 antibody



Immunocytochemistry analysis of SNX5 (green) in HeLa using SNX5 antibody, and DAPI (blue).

注意事项 (Note)

产品名称: SNX5 Rabbit Monoclonal Antibody
产品货号: AMRe85158



For research use only .