

产品名称: Rab5C Rabbit Monoclonal Antibody
产品货号: AMRe85163



产品概述 (Summary)

产品名称 (Production Name)	Rab5C Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide,0.05%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Rab5C
别名 (Alternative Names)	Ras-related protein Rab-5C; RABL; L1880; RAB5L; RAB5CL
基因 ID (Gene ID)	5878.0
蛋白 ID (SwissProt ID)	P51148.Recombinant protein of human RAB5C

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)


Protein transport.

产品名称: Rab5C Rabbit Monoclonal Antibody
产品货号: AMRe85163




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of RAB5C in K562, C6, 3T3, HeLa lysates using RAB5C antibody.



Immunocytochemistry analysis of Rab5C (green) in A549 cells using Rab5C antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using RAB5C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Rab5C Rabbit Monoclonal Antibody
产品货号: AMRe85163



注意事项 (Note)

For research use only .