

产品名称: TREX1 Rabbit Monoclonal Antibody
产品货号: AMRe85166



产品概述 (Summary)

产品名称 (Production Name)	TREX1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TREX1
别名 (Alternative Names)	AGS1; AGS5; CRV; DNase III; DRN3; HERNS; TREX1
基因 ID (Gene ID)	11277.0
蛋白 ID (SwissProt ID)	Q9NSU2.Recombinant protein of human TREX1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 33 kDa

研究背景 (Background)

TREX1 is the major 3'->5' DNA exonuclease in human cells. The protein is a non processive exonuclease that may serve a


产品名称: TREX1 Rabbit Monoclonal Antibody
产品货号: AMRe85166




proofreading function for a human DNA polymerase.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of TREX1 in HeLa lysates using TREX1 antibody.



Immunocytochemistry analysis of TREX1 (green) in HeLa using TREX1 antibody, and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using TREX1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: TREX1 Rabbit Monoclonal Antibody
产品货号: AMRe85166



注意事项 (Note)

For research use only .