

产品名称: Transferrin Receptor 1 Rabbit Monoclonal Antibody
产品货号: AMRe85181



产品概述 (Summary)

产品名称 (Production Name)	Transferrin Receptor 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Transferrin Receptor 1
别名 (Alternative Names)	TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD71
基因 ID (Gene ID)	7037.0
蛋白 ID (SwissProt ID)	P02786.A synthetic peptide of human Transferrin Receptor

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 85 kDa; Observed MW: 90 kDa

研究背景 (Background)


产品名称: Transferrin Receptor 1 Rabbit Monoclonal Antibody
产品货号: AMRe85181



Binds the telomeric double-stranded TTAGGG repeat and negatively regulates telomere length. Involved in the regulation of the mitotic spindle. Component of the shelterin complex (telosome) that is involved in the regulation of telomere length and protection.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Transferrin Receptor in HeLa lysates using Transferrin Receptor 1 antibody.

注意事项 (Note)

For research use only .