

产品名称: PMS2 Rabbit Monoclonal Antibody
产品货号: AMRe85182



产品概述 (Summary)

产品名称 (Production Name)	PMS2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PMS2
别名 (Alternative Names)	DNA mismatch repair gene; DNA mismatch repair protein PMS2; HNPCC4; PMS1 protein homolog 2
基因 ID (Gene ID)	5395.0
蛋白 ID (SwissProt ID)	P54278.A synthetic peptide of human PMS2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 96 kDa; Observed MW: 96 kDa

研究背景 (Background)


产品名称: PMS2 Rabbit Monoclonal Antibody
产品货号: AMRe85182



Assembly of the MutL-MutS-heteroduplex ternary complex in presence of RFC and PCNA is sufficient to activate endonuclease activity of PMS2. It introduces single-strand breaks near the mismatch and thus generates new entry points for the exonuclease EXO1 to degrade the strand containing the mismatch. DNA methylation would prevent cleavage and therefore assure that only the newly mutated DNA strand is going to be corrected.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PMS2 in Jurkat, HeLa lysates using PMS2 antibody.

注意事项 (Note)

For research use only .