

产品名称: StAR Rabbit Monoclonal Antibody
产品货号: AMRe85190



产品概述 (Summary)

产品名称 (Production Name)	StAR Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	StAR
别名 (Alternative Names)	STARD1
基因 ID (Gene ID)	6770.0
蛋白 ID (SwissProt ID)	P49675.A synthetic peptide of human StAR

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 32 kDa; Observed MW: 32 kDa

研究背景 (Background)

Steroidogenic acute regulatory protein (StAR) plays a significant role in cholesterol transport from the cytoplasmic outer


产品名称: StAR Rabbit Monoclonal Antibody
产品货号: AMRe85190



membrane to the inner mitochondrial membrane.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of StAR in K562, 3T3, HeLa lysates using StAR antibody.

注意事项 (Note)

For research use only .